



BIOLOGICAL INDICATOR INCUBATOR BIMA-111

BIOLOGICAL INDICATOR INCUBATOR BIMA-111

BIOLOGICAL INDICATOR INCUBATOR

Biological indicator incubators are designed for temperature incubation of the EO biological indicators. Both of them are compact, tabletop incubators.



LCD display.

Users can set the incubating temperature and time.

Exact correspondence of the set and actual temperature.

It is easy to set up and use.

The Incubator has many wells for placement of the biological indicator vials and the wells can be custom-made.

SPECIFICATIONS

Model	BIMA-111
Temp. control range	RT. +5 °C~80 °C
Timer range	1 min~99 h 59 min
Temp. control accuracy	±0.5 °C
Display accuracy	±0.1 °C
Heating time (25 to 57 °C)	≤6 min
Dimension (DxWxH) mm	110x150x80 mm
Weight	0.5 kg
Temp. stability @57 °C	±0.35 °C
Alt Name	Biological Indicator Incubator

ACCESSORIES FOR PURCHASE

No	Name	Description
1	Block	15xΦ10 mm
2	Block	15xΦ12 mm
3	Block	6xΦ15 mm
4	Block	6xΦ16 mm
5	Block	24xΦ8.8 mm
6	Block	18xΦ9 mm



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com