



## AIR-JACKETED MULTI GAS INCUBATOR BIMA-113

# AIR-JACKETED MULTI GAS INCUBATOR BIMA-113

## BIO-INCUBATOR

The Biological Indicator Incubator is compact, versatile, reliable and easy to use. Designed for single temperature incubation of the EO Biological Indicators. The thermometer is included with the incubator, providing easy verification of incubation temperature.



Compact, fits in the palm of your hand  
LCD display .it s easy to setup and use  
Different kinds of aluminum blocks are available  
Clear cover ensures temperature uniformity  
Portable;Optional 24V adapter

## SPECIFICATIONS

Model	BIMA-113
Operating temperature range	From 5°C above ambient to 100°C
Time range	1 min - 99 h 59 min
Temp. control accuracy	$\leq \pm 0.5^{\circ}\text{C}$
Display accuracy	0.1°C
Heating time	< 6 min (from 20°C to 57°C)
Temp. control accuracy @57°C	$\leq \pm 0.35^{\circ}\text{C}$
Dimensions (WxDxH, mm)	120 x 152 x 112
MOQ-5PCS	
Alt Name	Bio-Incubator

## ACCESSORIES FOR PURCHASE

No	Name	Specification
1	tube block	Ø10 mm x 32 mm x 15 (MOQ-5PCS)
2	tube block	Ø12 mm x 32 mm x 12 (MOQ-5PCS)
3	tube block	Ø15 mm x 32 mm x 8 (MOQ-5PCS)
4	tube block	Ø16 mm x 32 mm x 8 (MOQ-5PCS)
5	tube block	Ø9 mm x 18 (MOQ-5PCS)



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)