



PH METER BJB1B1 (BMET-201)

PH METER BJG1B1

BENCHTOP PH METER

The pH scale typically ranges from 0 to 14. According to the Nernst equation pH values are affected by temperature, hence the temperature compensation is a vital function of any pH Meter to ensure accurate pH measurement.



SPECIFICATIONS

Model	BJG1B1
Old Model	BMET-201
pH	
- Range	-2.00 to 14.00 pH
- Resolution	0.01 pH
- Accuracy	± 0.01 pH
- Calibration Points	Up to 2
- Standard Recognition	NIST buffers
mV	
Conductivity	
- Reference Temperature	/
Measurement	
- Reading Mode	Continuous
- Reading Prompts	Reading, Stable
- Temp. Compensation	MTC
Inputs	
- pH Electrode	BNC(Q9)
- Conductivity Probe	/
Display Options	
- Backlight	Yes
- Auto Shutdown	300, 600, 1200, 1800, 3600 sec., off
- IP Rating	IP54
General	
- Power	AC Adapter, 100-240 V AC input, DC 9V output
- Dimensions	200x160x63 mm
- Weight	600 g (1.32 lb)
Alt Name	Benchtop pH Meter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com