



MULTI-PARAMETER ANALYZER BJI1F1 (BMET-1105)

MULTI-PARAMETER ANALYZER BJ11F1

MULTIMETER

Multi-parameter Analyzer simultaneously measures Ion, Conductivity, Dissolved Oxygen, Temperature that otherwise requires four different meters. With the appropriate probes installed, it can also simultaneously measure several or one parameters.



SPECIFICATIONS

Model	BJ11F1
Old Model	BMET-1105
DC Voltage	
- 200 mV / 2 V / 20 V / 200 V	$\pm (0.5\% + 3)$
- 600 V	$\pm (0.8\% + 5)$
AC Voltage	
- 200 V / 600 V	$\pm (0.8\% + 3)$
DC Current	
- 200 μ A / 2000 μ A / 20 mA / 200 mA	$\pm (1.8\% + 5)$
AC Current	
Resistance	
- 200 Ω / 2 k Ω / 20 k Ω / 200 k Ω / 2 M Ω / 10 M Ω	$\pm (1.0\% + 3)$
- 20 M Ω	$\pm (1.2\% + 15)$
Display	2000 Counts
Data Hold	Yes
Continuity	Yes
NCV (Non-Contact Voltage)	Yes
Work Light	Yes
Power Supply	2 x 1.5 V AAA Batteries
Weight	Approx. 133.5 g / 146.5 g
Size	130 x 62 x 27 mm
Safety Standard / Rating	EN61010-1, EN61010-2-033, EN61326, CAT.III 600 V
Alt Name	Multimeter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com