



## MULTI-PARAMETER ANALYZER BJI1G1 (BMET-1106)

# MULTI-PARAMETER ANALYZER BJ11G1

## MULTIMETER

Multi-parameter Analyzer simultaneously measures Ion, Conductivity, Dissolved Oxygen, Temperature that otherwise requires four different meters. With the appropriate probes installed, it can also simultaneously measure several or one parameters.



## SPECIFICATIONS

| Model                                                          | BJ11G1                                           |
|----------------------------------------------------------------|--------------------------------------------------|
| Old Model                                                      | BMET-1106                                        |
| DC Voltage                                                     |                                                  |
| - 200 mV / 2 V / 20 V / 200 V                                  | $\pm (0.5\% + 2)$                                |
| - 600 V                                                        | $\pm (0.8\% + 2)$                                |
| AC Voltage                                                     |                                                  |
| - 2 V / 20 V / 200 V                                           | $\pm (1.2\% + 10)$                               |
| DC Current                                                     | $\pm (1.0\% + 2)$                                |
| - 200 $\mu$ A / 2 mA / 20 mA                                   | $\pm (1.0\% + 2)$                                |
| - 200 mA                                                       | $\pm (1.5\% + 2)$                                |
| - 10 A                                                         | $\pm (3.0\% + 2)$                                |
| Resistance                                                     |                                                  |
| - 200 $\Omega$ / 2 k $\Omega$ / 20 k $\Omega$ / 200 k $\Omega$ | $\pm (0.8\% + 2)$                                |
| - 2 M $\Omega$                                                 | $\pm (1.0\% + 2)$                                |
| Display                                                        | 2000 Counts                                      |
| Data Hold                                                      | Yes                                              |
| Diode Test                                                     | Yes                                              |
| Continuity                                                     | Yes                                              |
| Power Supply                                                   | 1 x 9V (6F22 Battery)                            |
| Weight                                                         | Approx. 169.5 g / 173 g / 174.5 g                |
| Size                                                           | 148 x 74 x 42 mm                                 |
| Safety Standard / Rating                                       | EN61010-1, EN61010-2-033, EN61326, CAT.III 600 V |
| Alt Name                                                       | Multimeter                                       |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)