



MULTI-PARAMETER ANALYZER BJI1H1 (BMET-1107)

MULTI-PARAMETER ANALYZER BJ11H1

MULTIMETER

Multi-parameter Analyzer simultaneously measures Ion, Conductivity, Dissolved Oxygen, Temperature that otherwise requires four different meters. With the appropriate probes installed, it can also simultaneously measure several or one parameters.



SPECIFICATIONS

Model	BJ11H1
Old Model	BMET-1107
DC Voltage	
- 200 mV / 2 V / 20 V / 200 V	$\pm (0.5\% + 2)$
- 600 V	$\pm (0.8\% + 2)$
AC Voltage	
- 2 V / 20 V / 200 V	$\pm (1.0\% + 3)$
DC Current	
- 200 μ A / 2 mA / 20 mA / 200 mA	$\pm (1.5\% + 3)$
AC Current	
DC (10 A Range)	
- 10 A	$\pm (3.0\% + 3)$
AC (10 A Range)	
Resistance	
- 200 Ω / 2 k Ω / 20 k Ω / 200 k Ω / 2 M Ω	$\pm (0.8\% + 2)$
- 20 M Ω	$\pm (1.0\% + 2)$
Temperature	
- -20°C to 1000°C	$\pm (3.0\% + 3)$
- 0°F to 1800°F	$\pm (3.0\% + 3)$
Display	2000 Counts
Data Hold	Yes
Continuity	Yes
Diode Test	Yes
Maximum Value Hold	Yes
Power Supply	1 x 9V (6F22 Battery)
Weight	Approx. 172.5 g / 173 g
Size	148 x 74 x 42 mm
Safety Standard / Rating	EN61010-1, EN61010-2-033, EN61326, CAT III 600 V
Alt Name	Multimeter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com