



## MULTI-PARAMETER ANALYZER BJI1J1 (BMET-1109)

# MULTI-PARAMETER ANALYZER BJI1J1

## MULTIMETER

Multi-parameter Analyzer simultaneously measures Ion, Conductivity, Dissolved Oxygen, Temperature that otherwise requires four different meters. With the appropriate probes installed, it can also simultaneously measure several or one parameters.



## SPECIFICATIONS

Model	BJI1J1
Old Model	BMET-1109
DC Voltage	
- 6 V / 60 V / 600 V	± (0.8% + 3)
AC Voltage	
- 60 mV / 600 mV / 6 V / 60 V / 600 V / 750 V	± (0.8% + 3)
Resistance	
- 6 kΩ / 60 kΩ / 600 kΩ	± (1.0% + 3)
- 6 MΩ / 10 MΩ	± (1.2% + 5)
Capacitance	
60 nF / 600 nF / 6 μF / 60 μF / 600 μF / 6 mF	± (4.0% + 5)
Frequency	
30 ~ 1000 Hz	± (0.5% + 2)
Display	6000 Counts
Auto Range	Yes
Bargraph Display	Yes
Continuity	Yes
NCV (Non-Contact Voltage)	Yes
Live Test	Yes
Data Hold	Yes
Work Light	Yes
Back Light	Yes
Auto Power Off	Yes
Power Supply	3 x 1.5 V AAA Batteries
Weight	Approx. 243 g
Size	169 x 83 x 53 mm
Safety Standard / Rating	EN61010-1, EN61010-2-033, EN61326, CAT.III 600 V
Alt Name	Multimeter



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)