



MULTI-PARAMETER ANALYZER BJI1M1 (BMET-1112)

MULTI-PARAMETER ANALYZER BJI1M1

MULTIMETER

Multi-parameter Analyzer simultaneously measures Ion, Conductivity, Dissolved Oxygen, Temperature that otherwise requires four different meters. With the appropriate probes installed, it can also simultaneously measure several or one parameters.



SPECIFICATIONS

Model	BJI1M1
Old Model	BMET-1112
DC Voltage	
- 400 mV / 4 V / 40 V / 400 V	$\pm (0.8\% + 2)$
- 600 V	$\pm (0.1\% + 2)$
AC Voltage	
- 40 V / 400 V / 600 V	$\pm (1.2\% + 10)$
DC Current	
- 400 μ A / 4 mA / 40 mA / 400 mA	$\pm (0.1\% + 2)$
- 400 mA	$\pm (1.5\% + 2)$
- 10 A	$\pm (3.0\% + 2)$
Resistance	
- 400 Ω / 4 k Ω / 40 k Ω / 400 k Ω	$\pm (0.8\% + 3)$
Display	4000 Counts
Manual Range	Yes
Bargraph Display	Yes
Data Hold	Yes
Diode Test	Yes
Continuity	Yes
NCV (Non-Contact Voltage Detection)	Yes
Back Light	Yes
Flashlight	Yes
Power Supply	2 x 1.5 V AAA
Weight	268.5 g
Size	166 x 78 x 48 mm
Safety Standard / Rating	EN61010-1, EN61010-2-033, EN61326, CAT.III 600 V
Alt Name	Multimeter



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com