



SLIDE CABINET BJJ1F1

SLIDE CABINET BJJ1F1

SLIDE DRYING CABINET



SPECIFICATIONS

Model	BJJ1F1
Base Material	SPCC 1.2 mm cold-roll steel sheets
Main Body and Drawers	SPCC 0.8 mm cold-roll steel sheets
Configuration	Base with three modules on top containing 18 drawers (6 drawers per module)
Capacity	Accommodates up to 9,000 microscope slides
Drawer Features	Numbered labels inside drawers make filing more convenient
Handles	Hidden handles made of ABS
Slot Manufacturing	Specifically manufactured using a metal mold
Labeling Slots	Manufactured through integrated stamping and forming operation
Dimensions (WxDxH)	450x478x1295 mm
Alt Name	Slide Drying Cabinet



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com