



## WHEELCHAIR BJP1B1 (BHBD-1201)

# WHEELCHAIR BJP1B1

## ELECTRIC TOILET LIFT

Electric Toilet lift is mainly suitable for the elderly, pregnant women, disabled and other people who have inconvenient movements



Adjustable Lifting: Adjustable lifting angle for seat height comfort at different heights

Strong Thrust: Double push strength up to 200kg, satisfying 150kg weight and sitting freely

Safe Voltage :DC 24V safe voltage, ultra-low energy consumption

Silent Operation: The whole process of silent operation ensures the comfort of sitting

Self-protection System: When a unilateral failure occurs, the machine automatically starts the self-protection function and stops running

One-button Reset: After use, long press the switch key to reset with one button, while locking other buttons to ensure safe use

## SPECIFICATIONS

|                   |                                                                                |
|-------------------|--------------------------------------------------------------------------------|
| Model             | BJP1B1                                                                         |
| Old Model         | BHBD-1201                                                                      |
| Material          | Iron, Plastic                                                                  |
| Seat height       | From 45 to 75 cm                                                               |
| Weight Capacity   | 150 kg                                                                         |
| Noise             | Almost silent electric motor when operation                                    |
| Inclination angle | 15°-16°                                                                        |
| Size              | 57 cm width, 65 cm length, 47 cm height                                        |
| For person        | Bariatric person, elderly person, patients disabled person, and pregnant women |
| Rated power       | 96 W / 2 A                                                                     |
| Voltage           | DC 24 V                                                                        |
| Power source      | Electric power                                                                 |
| Alt Name          | Electric Toilet lift                                                           |

## APPLICATIONS

Nursing home, rehabilitation center, Hospital



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)