



WHEELCHAIR BJP1B3 (BHBD-1203)

WHEELCHAIR BJP1B3

ELECTRIC TOILET LIFT

Electric Toilet lift is mainly suitable for the elderly, pregnant women, disabled and other people who have inconvenient movements



Adjustable Lifting: Adjustable lifting angle for seat height comfort at different heights

Strong Thrust: Double push strength up to 200kg, satisfying 150kg weight and sitting freely

Safe Voltage :DC 24V safe voltage, ultra-low energy consumption

Silent Operation: The whole process of silent operation ensures the comfort of sitting

Self-protection System: When a unilateral failure occurs, the machine automatically starts the self-protection function and stops running

One-button Reset: After use, long press the switch key to reset with one button, while locking other buttons to ensure safe use

SPECIFICATIONS

Model	BJP1B3
Old Model	BHBD-1203
Material	Iron, Plastic
Seat height	From 45 to 75 cm
Weight Capacity	150 kg
Noise	Almost silent electric motor when operation
Inclination angle	15°-16°
Size	57 cm width, 65 cm length, 47 cm height
For person	Bariatric person, elderly person, patients disabled person, and pregnant women
Rated power	96 W / 2 A
Voltage	DC 24 V
Power source	Electric power
Alt Name	Electric Toilet lift

APPLICATIONS

Nursing home, rehabilitation center, Hospital



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com