



## WHEELCHAIR BJP1C1 (BHBD-1101)

# WHEELCHAIR BJP1C1

## ELECTRIC LIFT PATIENT TRANSFER WHEELCHAIR

A wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.



Thickened aluminium alloy body  
strong bearing capacity  
180° opening and closing, easy to realize the transfer of various scenes  
Electric motor, one-button control, easy lifting operation, convenient and quick displacement  
Soft and comfortable armrest, comfortable to use  
360-degree universal wheel, with locking mechanism, shifting is more secure  
PU leather seat cushion, waterproof, can be showered, hollow design, a toilet can be placed below, Easy to extract

## SPECIFICATIONS

Model	BJP1C1
Old Model	BHBD-1101
Material	Iron, Plastic
Maximum loading weight	150 kg
Lifting Range	25 cm, from 40 cm to 65 cm.
Wheel	Two front wheels are with brake
Door width, chair can pass it	At least 55 cm
It suites for bed	Height of bed from 11 cm to 63 cm
Waterproof level	IP44
Dimensions	123x72.5x54.5 cm
Rated power	96 W
Voltage	DC 24 V
Power supply	Battery, rechargeable
Alt Name	Electric lift patient transfer Wheelchair

## APPLICATIONS

Patient transfer, patient lift, toilet-bath chair, Patient Indoor-outdoor locomotion



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)