



BACTERIAL COLONY COUNTER BJS1B1 (BCLC-103)

BACTERIAL COLONY COUNTER BJS1B1

MICROBIAL COLONY SPREADER



1. Four speeds options for user's requirement: 200RPM, 250RPM, 300RPM, 350RPM
2. Aluminum alloy disc can be sterilized by disinfectants and UV;
3. Continuous operation is for daily use;
4. Feasible 70/90mm plate, no need to replace the adapter

SPECIFICATIONS

Model	BJS1B1
Old Model	BCLC-103
Main power supply voltage	Voltage: ~220 V; Frequency: 50 Hz
Current	0.5 A
Diameter	170 mm (6.7 inches)
Height	80 mm (3.15 inches)
Weight	1.2 kg (2.6 lbs.)
Running Temperature	15 °C to 30 °C (59 °F to 95 °F)
Storage Temperature	-20 °C to 55 °C (-4 °F to 135 °F)
Relative Humidity	31 °C, ≤ 80%
Altitude	≤ 2000 m
Pollution Level	2
Voltage Classification	II
Alt Name	Microbial Colony Spreader

ITEMS INCLUDED

No	Name
1	Micro-Rot Alpha



1



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com