



BENCHTOP AUTOCLAVE BJT-601

BENCHTOP AUTOCLAVE BJT-601

TABLE TOP STEAM STERILIZER (CLASS N)



High-quality stainless steel inner container and sealing cover.
Microcomputer control with optional settings of 121°C and 134°C.
Touch buttons.
Equipped with a standard test interface.
Built-in steam generator for rapid supply of saturated steam.
Automatic protection functions: over-temperature protection, over-pressure protection, low water level protection, anti-dry heating.
Simple operation, safe and reliable.
European Class N standard.

SPECIFICATIONS

Model	BJT-601
Capacity	18L
Sterilization Temperature	134/121°C
Working Pressure	0.22 MPa
Features	
Power	1.5KW
Voltage	AC 110V 60HZ or AC 220V 50HZ
Design Pressure	0.25MPa
Design Temperature	139°C
Chamber Size (mm)	Φ250x360
Overall Dimensions (mm)	510x385x355
Packaging Dimensions (LxWxH) (mm)	590x470x410
Gross weight/Net weight	29/24kg
Alt Name	Table Top Steam Sterilizer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com