




LABORATORY VERTICAL AUTOCLAVE BJT1B2 (BAVT-2902)

LABORATORY VERTICAL AUTOCLAVE BJT1B2

VERTICAL AUTOCLAVE



Three times pre-vacuum
High quality stainless steel material
Microcomputer automatic control, arbitrary setting sterilization parameters
Safety interlocking device
LED touch screen display
Standard test interface
Self-expanding seal
With drying function
With stainless steel bucket or basket
With automatic protection function/over temperature protection; over pressure self discharge protection ; low water level protection; anti dry burning.
Automatic water

SPECIFICATIONS

Model	BJT1B2
Old Model	BAVT-2902
Capacity	50 L
Power	3.5 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C~134°C
Sterilization time setting range	4 ~ 120 min
Drying time setting range	0 ~ 240 min
Chamber volume	Ø 400x505
Bucket dimension	Ø 380x430
Basket dimension	Ø 360x280x1
Outer dimension	630x500x1180
Packing size	800x600x1230
Net weight	118 Kg
Gross Weight	138 Kg
Alt Name	Vertical Steam Sterilizer

FEATURES

Three times pre-vacuum
High quality stainless steel material
Microcomputer automatic control, arbitrary setting sterilization parameters
Safety interlocking device
LED touch screen display
Standard test interface
Self-expanding seal
With drying function
With stainless steel bucket or basket
With automatic protection function/over temperature protection; over pressure self discharge protection ; low water level protection;

anti dry burning.
Automatic water
External steam Generator
Buzzer reminder after sterilization, automatic stop.
Automatic discharge of cold air, automatic steam exhaust after sterilization
Internal printer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com