



LABORATORY VERTICAL AUTOCLAVE BJT1D2 (BAVT-3102)

LABORATORY VERTICAL AUTOCLAVE BJT1D2

VERTICAL AUTOCLAVE



Hand wheel translation and quick opening structure
High quality stainless steel material
Microcomputer automatic control, arbitrary setting sterilization parameters
Yr Safety interlocking device
Digital display and Touch buttons
Standard test interface
Self-expanding seal
With drying function
With automatic protection function: over temperature protection; over pressure protection; low water level protection; anti dry burning
With stainless steel bucket or basket Buzzer reminder after sterilization, automatic stop

SPECIFICATIONS

Model	BJT1D2
Old Model	BAVT-3102
Capacity	50 L
Power	3.5 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C~134°C
Sterilization time setting range	4 ~ 120 min
Drying time setting range	0 ~ 240 min
Chamber volume	Ø 350x525
Bucket dimension	Ø 330x460
Basket dimension	Ø 320x240x2
Outer dimension	480x480x1120
Packing size	580x630x1230
Net weight	80 Kg
Gross Weight	106 Kg
Alt Name	Vertical Steam Sterilizer

FEATURES

Hand wheel translation and quick opening structure
High quality stainless steel material
Microcomputer automatic control, arbitrary setting sterilization parameters
Yr Safety interlocking device
Digital display and Touch buttons
Standard test interface
Self-expanding seal
With drying function
With automatic protection function: over temperature protection; over pressure protection; low water level protection; anti dry burning

With stainless steel bucket or basket Buzzer reminder after sterilization, automatic stop
Automatic discharge of cold air,automatic steam exhaust after sterilization



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com