



LABORATORY VERTICAL AUTOCLAVE BJT1D3 (BAVT-3103)

LABORATORY VERTICAL AUTOCLAVE BJT1D3

VERTICAL AUTOCLAVE



Hand wheel translation and quick opening structure
High quality stainless steel material
Microcomputer automatic control, arbitrary setting sterilization parameters
Yr Safety interlocking device
Digital display and Touch buttons
Standard test interface
Self-expanding seal
With drying function
With automatic protection function: over temperature protection; over pressure protection; low water level protection; anti dry burning
With stainless steel bucket or basket Buzzer reminder after sterilization, automatic stop

SPECIFICATIONS

Model	BJT1D3
Old Model	BAVT-3103
Capacity	75 L
Power	4.5 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C~134°C
Sterilization time setting range	4 ~ 120 min
Drying time setting range	0 ~ 240 min
Chamber volume	Ø 400x625
Bucket dimension	Ø 380x560
Basket dimension	Ø 360x280x2
Outer dimension	525x525x1190
Packing size	630x680x1290
Net weight	90 Kg
Gross Weight	120 Kg
Alt Name	Vertical Steam Sterilizer

FEATURES

Hand wheel translation and quick opening structure
High quality stainless steel material
Microcomputer automatic control, arbitrary setting sterilization parameters
Yr Safety interlocking device
Digital display and Touch buttons
Standard test interface
Self-expanding seal
With drying function
With automatic protection function: over temperature protection; over pressure protection; low water level protection; anti dry burning

With stainless steel bucket or basket Buzzer reminder after sterilization, automatic stop
Automatic discharge of cold air,automatic steam exhaust after sterilization



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com