



LABORATORY VERTICAL AUTOCLAVE BJT1E6 (BAVT-3306)

LABORATORY VERTICAL AUTOCLAVE BJT1E6

VERTICAL AUTOCLAVE



Bolt structure
High quality stainless steel material
Microcomputer automatic control, arbitrary set of sterilization parameters
Digital display and touch buttons
Standard test interface
The container with automatic protection function: over temperature protection;
over pressure protection self discharge; low water level protection, anti dry burning
With stainless steel bucket or basket
Buzzer reminder after sterilization, automatic stop.
Automatic cooling air, sterilization end automatic exhaust steam

SPECIFICATIONS

Model	BJT1E6
Old Model	BAVT-3306
Capacity	150 L
Power	6.0 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C~134°C
Sterilization time setting range	4 ~ 120 min
Chamber volume	Ø 500x780
Bucket dimension	Ø 470x720
Basket dimension	Ø 460x350x2
Outer dimension	605x649x1210
Packing size	750x700x1370
Net weight	90 Kg
Gross Weight	116 Kg
Alt Name	Vertical Steam Sterilizer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com