



CLASS B AUTOCLAVE BJU-10203

CLASS B AUTOCLAVE BJU-10203

CLASS B BENCHTOP AUTOCLAVE



Complies with European CLASS B standard for effective sterilization.
Automatic door switch for convenient and safe operation.
Microcomputer control with LED display ensures precise and reliable sterilization.
Equipped with BD test, HELIX test, and Vacuum test for comprehensive sterilization performance.
Standard test interface for convenient testing and validation.
Automatic protection function includes over temperature, over pressure, low water level, anti dry burning, and safety door lock.
Safety interlocking device for secure operation.
Three times pre-vacuum enhances sterilization efficiency.
Buzzer reminder and automatic stop after sterilization cycle completion.

SPECIFICATIONS

Model	BJU-10203
Capacity	45L
Sterilization Temperature	134°C
Working Pressure	0.22 MPa
Features	
Consumption	3.5 kW
Power Supply	AC 220/110 V, 50/60 Hz
Design Pressure	0.25 MPa
Design Temperature	139°C
Rated Working Pressure	0.22 MPa
Rated Working Temperature	134°C
Water Supply	Built-in water tank
Chamber Size	Ø320 x 615 mm
External Size (WxDxH)	850 x 660 x 570 mm
Shipping Size (WxDxH)	980 x 770 x 720 mm
N.W./G.W.	127/102 kg
Alt Name	Class B Benchtop Autoclave

FEATURES

Complies with European CLASS B standard for effective sterilization.
Automatic door switch for convenient and safe operation.
Microcomputer control with LED display ensures precise and reliable sterilization.
Equipped with BD test, HELIX test, and Vacuum test for comprehensive sterilization performance.
Standard test interface for convenient testing and validation.
Automatic protection function includes over temperature, over pressure, low water level, anti dry burning, and safety door lock.
Safety interlocking device for secure operation.
Three times pre-vacuum enhances sterilization efficiency.
Buzzer reminder and automatic stop after sterilization cycle completion.
Automatic discharge of cold air and steam exhaust for efficient and safe operation.
Built-in steam generators provide quick production of saturated steam.
Standard USB Port and Printer for record-keeping and documentation.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com