



PATHOLOGY WORKSTATION BJU-3803

PATHOLOGY WORKSTATION BJU-3803

PATHOLOGY WORKSTATION

Workstations, accessories, and replacement parts for use during dissections, autopsies, necropsies, and other pathological procedures. Products include grossing/trimming stations and dissection tables, built-in sink and ventilation system.



SPECIFICATIONS

Model	BJU-3803
Type	Open type, no front glass door
External Size(WxDxH)(mm)	1800x760x1950
Working Surface Size(WxDxH)(mm)	1695x430x50
Work Surface Height	900mm
Exhaust Airflow Volume	1000-1200m ³ /h
Main Body	304 stainless steel
Work Surface	304 stainless steel
LED Lamp	36Wx1
UV Lamp	20Wx1
Standard Accessory	LED lampx1UV lampx1Spotlightx1Water-proof socketx2Adjustable hot and cold water tap and water sink(need to choose optional electrical water heater, then use hot water)Telescopic faucetStainless steel water sinkMagnetic tool holderActivated carbon filter4m exhaust PVC duct (Diameter: 300mm)
Optional Accessory	316L stainless steel work surface, thickness: 2 mmPulverizerCamera25x magnifierElectrical water heaterExternal exhaust blower
Power Supply	AC 220±10%V, 50/60Hz; 110±10%V, 60Hz
Consumption	700W
Shipping Size(WxDxH)(mm)	1980x910x2150
G.W.(kg)	380
Alt Name	Pathology Workstation



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com