



PASS BOX BJU-4501

PASS BOX BJU-4501

PASS BOX

Pass Box is used to transfer goods between rooms of different clean levels of clean room. It can reduce times of opening door, prevent the air convection between different rooms and minimize the pollution extent.



SPECIFICATIONS

Model	BJU-4501
External Size(WxDxH)(mm)	640x560x1000
Internal Size(WxDxH)(mm)	500x500x500
Material	304 Stainless steel
Interlock	Electronical
UV Lamp	8Wx1
HEPA Filter	Efficiency≥ 99.999% at 0.3μm
Clean Air Speed	0.36~0.54m/s
Gross Weight(kg)	75
Shipping Size(WxDxH)(mm)	700x700x1250
Alt Name	Air Shower Pass Box



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com