



PASS BOX BJU-4502

## PASS BOX BJU-4502

### PASS BOX

Pass Box is used to transfer goods between rooms of different clean levels of clean room. It can reduce times of opening door, prevent the air convection between different rooms and minimize the pollution extent.



## SPECIFICATIONS

| Model                    | BJU-4502                     |
|--------------------------|------------------------------|
| External Size(WxDxH)(mm) | 740x660x1100                 |
| Internal Size(WxDxH)(mm) | 600x600x600                  |
| Material                 | 304 Stainless steel          |
| Interlock                | Electronical                 |
| UV Lamp                  | 8Wx1                         |
| HEPA Filter              | Efficiency≥ 99.999% at 0.3μm |
| Clean Air Speed          | 0.36~0.54m/s                 |
| Gross Weight(kg)         | 95                           |
| Shipping Size(WxDxH)(mm) | 800x800x1350                 |
| Alt Name                 | Air Shower Pass Box          |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)