



PASS BOX BJU-4505

PASS BOX BJU-4505

PASS BOX

Pass Box is used to transfer goods between rooms of different clean levels of clean room. It can reduce times of opening door, prevent the air convection between different rooms and minimize the pollution extent.



SPECIFICATIONS

Model	BJU-4505
External Size(WxDxH)(mm)	740x660x680
Internal Size(WxDxH)(mm)	600x600x600
Material	304 Stainless steel
Interlock	Electronical
UV Lamp	8Wx1
HEPA Filter	-
Clean Air Speed	-
Gross Weight(kg)	70
Shipping Size(WxDxH)(mm)	800x800x850
Alt Name	Pass Box



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com