



PASS BOX BJU-4506

PASS BOX BJU-4506

PASS BOX

Pass Box is used to transfer goods between rooms of different clean levels of clean room. It can reduce times of opening door, prevent the air convection between different rooms and minimize the pollution extent.



SPECIFICATIONS

Model	BJU-4506
External Size(WxDxH)(mm)	840x760x780
Internal Size(WxDxH)(mm)	700x700x700
Material	304 Stainless steel
Interlock	Electronical
UV Lamp	8Wx1
HEPA Filter	-
Clean Air Speed	-
Gross Weight(kg)	80
Shipping Size(WxDxH)(mm)	900x900x950
Alt Name	Pass Box



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com