



DISPENSING BOOTH BJU-4701

DISPENSING BOOTH BJU-4701

DISPENSING BOOTH (SAMPLING OR WEIGHING BOOTH)

The dispensing booth is a local purification equipment dedicated to places such as pharmaceuticals, microbiological research and scientific experiments. It provides a kind of vertical, unidirectional airflow that produces negative pressure in the working area, part

of the clean air circulates in the working area, part is discharged to the nearby area to prevent cross-contamination, to ensure the high cleanliness in the working area.



SPECIFICATIONS

Model	BJU-4701
External Size(WxDxH)	1000x1350x2400mm
Internal Size (WxDxH)	900x850x1950mm
Air Velocity	0.38-0.58m/s
Supply air volume	1800m³/h
Exhaust air volume	10%-15%
HEPA filter size (WxDxH)	720x590x90mm
Clean Level	ISO 5 (Class 100)
Illumination	≥300 lux
Sound Emission	≤75dB(A)
Consumption	0.35KW
Power Supply	AC three-phase 380V/50Hz(Optional: AC single-phase 220/110V,50Hz)
N.W./G.W. (kg)	290/320
Shipping Size (WxDxH)(mm)	1200x1550x2700
Shipping Volume (m³)	5.1
Alt Name	Dispensing Booth (Sampling or Weighing Booth)



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com