



DISPENSING BOOTH BJU-4702

DISPENSING BOOTH BJU-4702

DISPENSING BOOTH (SAMPLING OR WEIGHING BOOTH)

The dispensing booth is a local purification equipment dedicated to places such as pharmaceuticals, microbiological research and scientific experiments. It provides a kind of vertical, unidirectional airflow that produces negative pressure in the working area, part of the clean air circulates in the working area, part is discharged to the nearby area to prevent cross-contamination, to ensure the high cleanliness in the working area.



SPECIFICATIONS

Model	BJU-4702
External Size(WxDxH)	1300x1350x2400mm
Internal Size (WxDxH)	1200x850x1950mm
Air Velocity	0.38-0.58m/s
Supply air volume	2600m³/h
Exhaust air volume	0%-15%
HEPA filter size (WxDxH)	1020x590x90mm
Clean Level	ISO 5 (Class 100)
Illumination	≥300 lux
Sound Emission	≤75dB(A)
Consumption	0.55KW
Power Supply	AC three-phase 380V/50Hz(Optional: AC single-phase 220/110V,50Hz)
N.W./G.W. (kg)	380/420
Shipping Size (WxDxH)(mm)	1500x1550x2700
Shipping Volume (m³)	6.3
Alt Name	Dispensing Booth (Sampling or Weighing Booth)



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com