



PRODUCT CATALOG



NITROGEN GENERATOR BJV-102

NITROGEN GENERATOR BJV-102

NITROGEN GENERATOR



Simple operation, just connect compressed air and turn on the power switch to generate high-purity nitrogen;
The electrolytic cell has a large electrolysis area, low cell temperature and long service life;
Equipped with a liquid-proof device to ensure that the instrument does not have liquid return;
Automatic emptying, power failure has the function of automatically closing the gas path;
Pressure and flow are automatically displayed, automatic constant pressure and constant flow, and the nitrogen flow can be fully automatically adjusted according to the amount used.

SPECIFICATIONS

Model	BJV-102
Output flow: ml/min	0-500
Input pressure: (MPa)	0.45
Output pressure: (MPa)	0.4
Gas purity: (%)	99.995
Input power: (W)	150
External Dimensions(W*D*H)(mm)	310*200*370
Package Dimension (W*D*H)(mm)	440*320*500
G.W.(kg)	11
Alt Name	Nitrogen Generator



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com