



PRODUCT CATALOG



HYDROGEN AIR GENERATOR BJV-201

HYDROGEN AIR GENERATOR BJV-201

HYDROGEN GAS GENERATOR



Simple operation, high-purity hydrogen can be produced by turning on the power supply;
High electrolysis efficiency, long-term use, purity and flow rate will not decay;
High reliability, simple daily use and maintenance, only need to add pure water regularly;
Pressure and flow rate are automatically displayed, automatic constant pressure and constant flow, and hydrogen flow rate can be fully automatically adjusted according to the amount used;
Equipped with safety devices, sensitive and reliable, and equipped with anti-return liquid devices to ensure that the instrument has no return liquid phenomenon.

SPECIFICATIONS

Model	BJV-201
Output flow: ml/min	0-300
Output pressure: (MPa)	0.3Mpa
Gas purity: (%)	99.999
Input power: (w)	150
External Dimensions(W*D*H)(mm)	310*200*370
Package Dimension (W*D*H)(mm)	430*310*480
G.W.(kg)	11.2
Alt Name	Hydrogen Gas Generator



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com