



# PRODUCT CATALOG



## BENCHTOP SHAKING INCUBATOR BJW-401

# BENCHTOP SHAKING INCUBATOR BJW-401

WRIST-ACTION FLASK SHAKER



Large backlit LCD display for clear parameter setting and real-time observation; password - protected interface prevents repeated or incorrect operations. High-quality servo motor providing precise speed control, excellent high-speed performance, and strong stability. Power-failure recovery to prevent data loss caused by unexpected outages or system errors. Streamlined exterior with electrostatic powder-coated finish for enhanced durability and corrosion resistance. Simulated wrist-action shaking delivers uniform and consistent oscillation. Multiple clamp options available to accommodate different flask sizes and application needs.

## SPECIFICATIONS

Model	BJW-401
No-load shaking speed	10-300rpm
Shaking speed accuracy	±1rpm
Shaking angle	18°
Display mode	Non-touch LCD screen
Timing range	0-999.9hours
Maximum capacity	250mlx8 or 500mlx8
Power supply	AC220±10%, 50 ~ 60Hz
Overall Dimensions(L*W*H)	792mmx332mmx310mm
Net weight	36.5kg
Noise level	Below 65dB
Alt Name	Wrist-Action Flask Shaker



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)