



## 2°C TO 8°C PHARMACY REFRIGERATOR BJX-1204

# 2°C TO 8°C PHARMACY REFRIGERATOR BJX-1204

## 2-8°C PHARMACY REFRIGERATOR



### Temperature Control System

Microprocessor control, Large LED display inner temperature clearly, and with easy viewing.

The inner temperature can be adjustable at the range of 2°C~8°C with an increment of 0.1.

Multiple sensors to display or control the temperature.

The Lowest and Highest temperature recorded value record.

### Safety Control

Malfunction alarms: high temperature alarm, low temperature alarm, Power failure alarm, Door open alarm.

Over temperature alarm system, set the alarm temperature as requirements.

Three alarming methods: buzzer alarming, digital flashing alarming and Remote al

## SPECIFICATIONS

|                            |                                |
|----------------------------|--------------------------------|
| Model                      | BJX-1204                       |
| Capacity                   | 255L                           |
| Temperature range          | 2~8°C                          |
| Cabinet Type               | Upright                        |
| Climate Class              | N                              |
| Cooling Type               | Forced air cooling             |
| Defrost Mode               | Auto                           |
| Refrigerant                | R600a / R290                   |
| Cooling performance(°C)    | 4                              |
| Controller                 | Microprocessor                 |
| Display                    | LED/---                        |
| Alarm                      | Audible, Remote                |
| Interior                   | Aluminum powder coating(White) |
| Exterior                   | PCM(White)                     |
| Power Supply(V/HZ)         | 220/50                         |
| Power(W)                   | 170                            |
| Power Consumption(kWh/24H) | 2.1                            |
| Net/Gross Weight(approx)   | 65kg/77kg                      |
| Interior Dimensions(WDH)   | 500x420x1290mm                 |
| Exterior Dimensions(WDH)   | 600x560x1525mm                 |
| Package Dimensions(WDH)    | 660x620x1660mm                 |
| High/Low Temperature       | Yes                            |
| Remote Alarm               | Yes                            |
| Sensor Failure             | Yes                            |
| Door Ajar                  | Yes                            |
| Lockage                    | Yes                            |
| Inner LED light            | Yes                            |
| Power Failure              | optional                       |
| Foot                       | Yes                            |
| Caster                     | None                           |
| Test Hole                  | Yes                            |

|                      |                             |
|----------------------|-----------------------------|
| Shelves              | 6                           |
| Foaming/Glass Door   | optional                    |
| USB Interface        | optional                    |
| Temperature Recorder | optional                    |
| Alt Name             | 2-8°C Pharmacy Refrigerator |

## FEATURES

### Temperature Control System

Microprocessor control, Large LED display inner temperature clearly, and with easy viewing.

The inner temperature can be adjustable at the range of 2°C~8°C with an increment of 0.1.

Multiple sensors to display or control the temperature.

The Lowest and Highest temperature recorded value record.

### Safety Control

Malfunction alarms: high temperature alarm, low temperature alarm, Power failure alarm, Door open alarm.

Over temperature alarm system, set the alarm temperature as requirements.

Three alarming methods: buzzer alarming, digital flashing alarming and Remote alarm.

Protection function: power-on delay, time between stops.

Temperature recorder monitor system (Optional).

### Refrigeration System

Highly efficient famous brand compressor and fan.

Forced-Air Circulation for large airflow with special air ducts to guarantee inner temperature consistency.

Large area of the evaporator design, cooling strong and high temperature uniformity, While effectively suppressing frost, reducing defrosting process.

CFC-Free/HCFC refrigerant.

### Ergonomic Design

Safety door lock, preventing unauthorized access.

Auto defrost function guarantees a good performance in high ambient temperature and humidity area.

Wide voltage design from 192V to 242V.

Inner LED lighting and caster design.

No condensation glass door.

Adjustable shelf design.

Foaming door is optional.

Stopper on rear against back wall.

## APPLICATIONS

Auto Defrost, Forced-Air Circulation Suitable for hospitals, blood banks, epidemic

prevention, animal husbandry areas, pharmaceutical companies and research institutes. Designed to store pharmaceuticals, medicine, vaccines, biological materials, testing reagents and laboratory materials.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)