



-86°C ULTRA LOW TEMPERATURE CHEST FREEZER BJX-303

-86°C ULTRA LOW TEMPERATURE FREEZER BJX-303

-86°C CHEST ULT FREEZER



SPECIFICATIONS

Model	BJX-303
Capacity	730L
Temperature range	-40~-86°C
Cabinet Type	Chest
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	R290/R170
Cooling performance(°C)	-80
Power (W)	980
Energy consumption (KW.H/24H)	14
Exterior Material	Galvanized steel powder coating
Interior Material	Stainless steel
Insulation Material	PUF+VIP
Net/Gross Weight(approx)	275kg/318kg
Interior Dimensions(WDH)	1310x740x730mm
Exterior Dimensions(WDH)	1965x1000x1035mm
Package Dimensions(WDH)	2040x1060x1190mm
Container Loading (20'/40')	10/20
Thickness of Cabinet Foamed Layer	90mm
Thickness of Door	80mm
Capacity for 2 inch boxes	540
Power Supply(V/Hz)	220/50
Display	7" Touch screen with USB port
High/Low Temperature	Yes
Hot Condenser	Yes
Sensor Failure	Yes
High Ambient Temp	Yes
Alarm mode	Sound and light alarm, remote alarm terminal
Caster	Yes
Test Hole	Yes
Temperature Recorder	Optional
Door locking device	Yes
Handle	Yes

Pressure balance hole	Yes
Racks & Boxes	Optional
Alt Name	-86°C Chest ULT Freezer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
 Email: info@biolabscientific.com | Website: www.biolabscientific.com