



-40°C CHEST FREEZER BJX-601

-40°C FREEZER BJX-601

-40°C CHEST DEEP FREEZER



SPECIFICATIONS

Model	BJX-601
Capacity	100L
Temperature range	-10~-40°C
Cooling performance(°C)	-40
Controller	Microprocessor
Display	LED
Interior	Aluminum powder coating
Exterior	Galvanized steel powder coating
Power Supply(V/HZ)	220/50
Power(W)	180
Refrigerant	R290
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Cabinet Type	Chest
Net/Gross Weight(approx)	40kg/47kg
Interior Dimensions(W*D*H)	520x345x615mm
Exterior Dimensions(W*D*H)	700x550x850mm
Package Dimensions(W*D*H)	750x595x940mm
High/Low Temperature	Yes
Sensor Error	Yes
Lockage	None
Caster	4
Foot	None
Test Hole	None
Baskets/Inner Doors	1/-
Temperature Recorder	Optional
Cryo racks	Optional
Alt Name	-40°C Chest Deep Freezer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com