



-40°C CHEST FREEZER BJX-605

-40°C FREEZER BJX-605

-40°C CHEST DEEP FREEZER



SPECIFICATIONS

Model	BJX-605
Capacity	470L
Temperature range	-10~-40°C
Cooling performance(°C)	-40
Controller	Microprocessor
Display	LED
Interior	Aluminum powder coating
Exterior	Galvanized steel powder coating
Power Supply(V/HZ)	220/50
Power(W)	350
Refrigerant	R404a
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Cabinet Type	Chest
Net/Gross Weight(approx)	110kg/130kg
Interior Dimensions(W*D*H)	1710x485x625mm
Exterior Dimensions(W*D*H)	1900x765x885mm
Package Dimensions(W*D*H)	1950x780x1060mm
High/Low Temperature	Yes
Sensor Error	Yes
Lockage	Yes
Caster	4
Foot	None
Test Hole	None
Baskets/Inner Doors	2/-
Temperature Recorder	Optional
Cryo racks	Optional
Alt Name	-40°C Chest Deep Freezer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com