



-25°C CHEST FREEZER BJB-903

-25°C FREEZER BJX-903

-25°C DEEP FREEZER



Safety Control

Malfunction alarms: high temperature alarm, low temperature alarm, Sensor failure.

Over temperature alarm system, set the alarm temperature as requirements.

Two alarming methods: buzzer alarming and digital flashing alarming.

Protection function: power-on delay, time between stops.

Refrigeration System

Highly efficient famous brand compressor and fan, with a high efficient refrigeration effect.

Optimized refrigeration technology, low noise, more energy efficient.

CFC-Free overall polyurethane foaming insulation.

CFC-Free refrigerant.

SPECIFICATIONS

Model	BJX-903
Capacity	285L
Temperature range	-10~-25°C
Cabinet Type	Chest
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	R290
Cooling performance(°C)	-25
Controller	Microprocessor
Display	LED
Interior	Aluminum powder coating
Exterior	Galvanized steel powder coating
Power Supply(V/HZ)	220/50
Power(W)	210
Net/Gross Weight(approx)	75kg/88kg
Interior Dimensions(W*D*H)	1080x485x625mm
Exterior Dimensions(W*D*H)	1270x765x885mm
Package Dimensions(W*D*H)	1310x780x1060mm
High/Low Temperature	Yes
Sensor failure	Yes
Lock	Yes
Caster	4
Foot	None
Test Hole	None
Baskets/Inner Doors	2/-
Temperature Recorder	Optional
Cryo racks	Optional
Alt Name	-25°C Deep Freezer

FEATURES

Safety Control

Malfunction alarms: high temperature alarm, low temperature alarm, Sensor failure.

Over temperature alarm system, set the alarm temperature as requirements.

Two alarming methods: buzzer alarming and digital flashing alarming.

Protection function: power-on delay, time between stops.

Refrigeration System

Highly efficient famous brand compressor and fan, with a high efficient refrigeration effect.

Optimized refrigeration technology, low noise, more energy efficient.

CFC-Free overall polyurethane foaming insulation.

CFC-Free refrigerant.

Ergonomic Design

Safety door lock, preventing unauthorized access.

Wide voltage design from 192V to 242V.

Pressure-balanced port design, easy to open the door.

Adjustable shelves design.

Heavy type caster and locking adjustable screw, easy to move and fix.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com