



## 2°C TO 8°C PHARMACY REFRIGERATOR BJY-1007

# 2°C TO 8°C PHARMACY REFRIGERATOR BJY-1007

MEDICAL REFRIGERATOR O SERIES



## Safety Design

Multiple alarm functions to protect your reagents in real time

## Real Visibility

Built-in LED lighting for easy observation

## Heating Glass Door

Ensure that the door body will not condense under high temperature and humidity environment

## Constant Temperature

Forces air cooling system keeps the temperature stable and even and keeps the refrigerator dry

## SPECIFICATIONS

|                   |                               |
|-------------------|-------------------------------|
| Model             | BJY-1007                      |
| Capacity          | 1360L                         |
| Temperature Range | 2~8°C                         |
| Outer size(W*D*H) | 1800*745*1870mm               |
| Inner size(W*D*H) | 1680*610*1380mm               |
| Alt Name          | Medical Refrigerator O series |

## APPLICATIONS

Used to store blood, reagents, medicines, etc.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)