



-60°C ULTRA LOW PORTABLE CAR FREEZER BJY-103

PORTABLE FREEZER BGY-103

PORTABLE CAR ULT FREEZER DC12V/DC24V (DC)



High Energy Core
Famous compressor, more powerful cooling
Anti-vibration & Anti-shaking
Anti-vibration and anti-shake design for different road conditions
Dual-use
Support both DC and AC working,maximize the protection you want

SPECIFICATIONS

Model	BJY-103
Capacity	20L
Temperature Range	-60°C
Outer size(W*D*H)	780*410*515mm
Inner size(W*D*H)	327*220*305mm
Alt Name	Portable Car ULT Freezer DC12V/DC24(DC)

APPLICATIONS

For preservation of viruses, bacteria, cells, skin tissues, biological materials,etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com