



2°C TO 8°C PHARMACY REFRIGERATOR BJY-1211

2°C TO 8°C PHARMACY REFRIGERATOR BJY-1211

MEDICAL REFRIGERATOR 0 SERIES



Safety Design

Multiple alarm functions to protect your reagents in real time.

Real Visibility

Built-in LED lighting for easy observation.

Heating Glass Door

Ensure that the door body will not condense under high temperature and humidity environment.

Constant Temperature

Forces air cooling system keeps the temperature stable and even and keeps the refrigerator dry.

SPECIFICATIONS

Model	BJY-1211
Capacity	1160L
Temperature Range	2~8°C
Outer size(W*D*H)	1800*745*1870mm
Inner size(W*D*H)	1680*610*1380mm
Alt Name	Medical Refrigerator 0 series

APPLICATIONS

Used to store blood, reagents, medicines, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com