



## 4°C BLOOD BANK REFRIGERATOR BJY-1407

# 4°C BLOOD BANK REFRIGERATOR BJY-1407

4±1°C MEDICAL REFRIGERATOR O SERIES



Safety Design  
Multiple alarm functions to protect your reagents in real time  
Real Visibility  
Built-in LED lighting for easy observation  
Heating Glass Door  
Ensure that the door body will not condense under high temperature and humidity environment  
Constant Temperature  
Forces air cooling system keeps the temperature stable and even and keeps the refrigerator dry

## SPECIFICATIONS

|                   |                                     |
|-------------------|-------------------------------------|
| Model             | BJY-1407                            |
| Capacity          | 660L                                |
| Temperature Range | 4±1°C                               |
| Outer size(W*D*H) | 1100*650*2080mm                     |
| Inner size(W*D*H) | 1112*560*1232mm                     |
| Alt Name          | 4±1°C Medical Refrigerator O series |

## APPLICATIONS

Used to store blood, reagents, medicines, etc.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)