



4°C BLOOD BANK REFRIGERATOR BJY-1503

4°C BLOOD BANK REFRIGERATOR BJY-1503

4±1°C MEDICAL REFRIGERATOR 8 SERIES



Safety Design

Multiple alarm functions to protect your reagents in real time

Real Visibility

Built-in LED lighting for easy observation

Heating Glass Door

Ensure that the door body will not condense under high temperature and humidity environment

Constant Temperature

Forces air cooling system keeps the temperature stable and even and keeps the refrigerator dry

SPECIFICATIONS

Model	BJY-1503
Capacity	668L
Temperature Range	4±1°C
Outer size(W*D*H)	1220*625*1885 mm
Inner size(W*D*H)	1100*504*1325 mm
Alt Name	4±1°C Medical Refrigerator 8 series

APPLICATIONS

Used to store blood, reagents, medicines, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com