



4°C BLOOD BANK REFRIGERATOR BJY-1504

4°C BLOOD BANK REFRIGERATOR BJY-1504

4±1°C MEDICAL REFRIGERATOR 8 SERIES



Safety Design
Multiple alarm functions to protect your reagents in real time
Real Visibility
Built-in LED lighting for easy observation
Heating Glass Door
Ensure that the door body will not condense under high temperature and humidity environment
Constant Temperature
Forces air cooling system keeps the temperature stable and even and keeps the refrigerator dry

SPECIFICATIONS

Model	BJY-1504
Capacity	868L
Temperature Range	4±1°C
Outer size(W*D*H)	1220*755*1885 mm
Inner size(W*D*H)	1100* 634*1325 mm
Alt Name	4±1°C Medical Refrigerator 8 series

APPLICATIONS

Used to store blood, reagents, medicines, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com