



-60°C ULTRA LOW TEMPERATURE UPRIGHT FREEZER BJY-1801

-60°C ULTRA LOW TEMPERATURE FREEZER BGY-1801

VERTICAL ULTRA LOW TEMPERATURE FREEZER



Efficiency & Reliability
Self-compounding system for superior cooling
Famous compressor with more stable performance
and longer using life
LED/LCD Display
Real-time display of operating status,
more intuitive and clear
Wide Voltage Design
Stable operation between 187V-242V. Freezer can be
worked stably in unstable voltage
Safety Design
Double lock, double protection

SPECIFICATIONS

Model	BGY-1801
Capacity	98L
Temperature Range	-60°C
Outer size(W*D*H)	620*620*1240mm
Inner size(W*D*H)	400*390*660mm
Alt Name	Vertical Ultra Low Temperature Freezer

APPLICATIONS

For preservation of viruses, bacteria, cells, skin tissues, biological materials, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com