



-60°C ULTRA LOW TEMPERATURE UPRIGHT FREEZER BJY-1802

-60°C ULTRA LOW TEMPERATURE FREEZER BJY-1802

VERTICAL ULTRA LOW TEMPERATURE FREEZER



Efficiency & Reliability
Self-compounding system for superior cooling
Famous compressor with more stable performance
and longer using life
LED/LCD Display
Real-time display of operating status,
more intuitive and clear
Wide Voltage Design
Stable operation between 187V-242V. Freezer can be
worked stably in unstable voltage
Safety Design
Double lock, double protection

SPECIFICATIONS

Model	BJY-1802
Capacity	188L
Temperature Range	-60°C
Outer size(W*D*H)	620*645*1572mm
Inner size(W*D*H)	440*425*985mm
Alt Name	Vertical Ultra Low Temperature Freezer

APPLICATIONS

For preservation of viruses, bacteria, cells, skin tissues, biological materials, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com