



## -60°C ULTRA LOW TEMPERATURE UPRIGHT FREEZER BJY-1804

# -60°C ULTRA LOW TEMPERATURE FREEZER BJY-1804

## VERTICAL ULTRA LOW TEMPERATURE FREEZER



Efficiency & Reliability  
Self-compounding system for superior cooling  
Famous compressor with more stable performance and longer using life  
LED/LCD Display  
Real-time display of operating status, more intuitive and clear  
Wide Voltage Design  
Stable operation between 187V-242V. Freezer can be worked stably in unstable voltage  
Safety Design  
Double lock, double protection

## SPECIFICATIONS

|                   |                                        |
|-------------------|----------------------------------------|
| Model             | BJY-1804                               |
| Capacity          | 368L                                   |
| Temperature Range | -60°C                                  |
| Outer size(W*D*H) | 690*705*1965mm                         |
| Inner size(W*D*H) | 510*515*1366mm                         |
| Alt Name          | Vertical Ultra Low Temperature Freezer |

## APPLICATIONS

For preservation of viruses, bacteria, cells, skin tissues, biological materials, etc.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)