



-60°C ULTRA LOW TEMPERATURE CHEST FREEZER BGY-2006

-60°C ULTRA LOW TEMPERATURE FREEZER BJY-2006

HORIZONTAL ULTRA LOW TEMPERATURE FREEZER 0 SERIES



Temperature Control

Adopting high-precision sensors, the temperature difference is controlled within $\pm 3^{\circ}\text{C}$

LED Display

Clear and intuitive real-time monitoring of operating status

Superior Performance

Famous compressor and fan, thickened foam layer make the refrigerator's performance better

Multiple alarm protections

Multiple alarm protections keep what's stored safer

SPECIFICATIONS

Model	BJY-2006
Capacity	520L
Temperature Range	-60°C
Outer size(W*D*H)	1855*890*1005mm
Inner size(W*D*H)	1220*600*695mm
Alt Name	Horizontal Ultra Low Temperature Freezer 0 series



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com