



# -80°C ULTRA LOW TEMPERATURE CHEST FREEZER BJY-2107

# -86°C ULTRA LOW TEMPERATURE FREEZER BJY-2107

## HORIZONTAL ULTRA LOW TEMPERATURE FREEZER 0 SERIES



### Temperature Control

Adopting high-precision sensors, the temperature difference is controlled within  $\pm 3^{\circ}\text{C}$

### LED Display

Clear and intuitive real-time monitoring of operating status

### Superior Performance

Famous compressor and fan, thickened foam layer make the refrigerator's performance better

### Multiple alarm protections

Multiple alarm protections keep what's stored safer

## SPECIFICATIONS

|                   |                                                   |
|-------------------|---------------------------------------------------|
| Model             | BJY-2107                                          |
| Capacity          | 620L                                              |
| Temperature Range | -80°C                                             |
| Outer size(W*D*H) | 2115*890*1005mm                                   |
| Inner size(W*D*H) | 1480*600*695mm                                    |
| Alt Name          | Horizontal Ultra Low Temperature Freezer 0 series |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)