



-105°C ULTRA LOW TEMPERATURE CHEST FREEZER BJY-2504

-105°C ULTRA LOW TEMPERATURE FREEZER BJY-2504

HORIZONTAL ULTRA LOW TEMPERATURE FREEZER 8 SERIES



Temperature Control

Adopting high-precision sensors, the temperature difference is controlled within $\pm 3^{\circ}\text{C}$

LED Display

Clear and intuitive real-time monitoring of operating status

Superior Performance

Famous compressor and fan, thickened foam layer make the refrigerator's performance better

Multiple alarm protections

Multiple alarm protections keep what's stored safer

SPECIFICATIONS

Model	BJY-2504
Capacity	358L
Temperature Range	-105°C
Outer size(W*D*H)	1810*895*1090mm
Inner size(W*D*H)	1080*500*660mm
Alt Name	Horizontal Ultra Low Temperature Freezer 8 series



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com