



## -25°C UPRIGHT FREEZER BGY-2806

## -25°C FREEZER BJY-2806

### VERTICAL ULTRA LOW TEMPERATURE FREEZER



Efficiency & Reliability  
Self-compounding system for superior cooling  
Famous compressor with more stable performance and longer using life  
LED/LCD Display  
Real-time display of operating status, more intuitive and clear  
Wide Voltage Design  
Stable operation between 187V-242V. Freezer can be worked stably in unstable voltage  
Safety Design  
Double lock, double protection

## SPECIFICATIONS

|                   |                                        |
|-------------------|----------------------------------------|
| Model             | BJY-2806                               |
| Capacity          | 308L                                   |
| Temperature Range | -25°C                                  |
| Outer size(W*D*H) | 810*760*1810                           |
| Inner size(W*D*H) | 570*485*1140                           |
| Alt Name          | Vertical Ultra Low Temperature Freezer |

## APPLICATIONS

For preservation of viruses, bacteria, cells, skin tissues, biological materials, etc.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)