



**-60°C ULTRA LOW PORTABLE CAR FREEZER
BGY-303**

PORTABLE FREEZER BJY-303

PORTABLE CAR ULT FREEZER (AC)



High Energy Core
Famous compressor, more powerful cooling
Anti-vibration & Anti-shaking
Anti-vibration and anti-shake design for different road conditions
Dual-use
Support both DC and AC working,maximize the protection you want

SPECIFICATIONS

Model	BJY-303
Capacity	20L
Temperature Range	-60°C
Outer size(W*D*H)	780*410*515mm
Inner size(W*D*H)	327*220*305mm
Alt Name	Portable Car ULT Freezer (AC)

APPLICATIONS

For preservation of viruses, bacteria, cells, skin tissues, biological materials,etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com