



-86°C ULTRA LOW TEMPERATURE UPRIGHT FREEZER BGY-3112

-86°C ULTRA LOW TEMPERATURE FREEZER BJY-3112

VERTICAL ULTRA LOW TEMPERATURE FREEZER



Efficiency & Reliability
Self-compounding system for superior cooling
Famous compressor with more stable performance and longer using life
LED/LCD Display
Real-time display of operating status, more intuitive and clear
Wide Voltage Design
Stable operation between 187V-242V. Freezer can be worked stably in unstable voltage
Safety Design
Double lock, double protection

SPECIFICATIONS

Model	BJY-3112
Capacity	938L
Temperature Range	-86°C
Outer size(W*D*H)	1240*1000*2000mm
Inner size(W*D*H)	1000*700*1335mm
Alt Name	Vertical Ultra Low Temperature Freezer

APPLICATIONS

For preservation of viruses, bacteria, cells, skin tissues, biological materials, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com