



2°C TO 8°C SMART VACCINE REFRIGERATOR BJY-3202

2°C TO 8°C PHARMACY REFRIGERATOR BJY-3202

SMART MEMORY VACCINE CABINET



SPECIFICATIONS

Model	BJY-3202
Capacity	360L
Temperature Range	2-8°C
Outer size(W*D*H)	680*660*1995mm
Inner size(W*D*H)	580*490*1400mm
Alt Name	Smart Memory Vaccine Cabinet



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com